



**FINISHED
ORGANIC RESIN LENSES
TECHNICAL DATA**

**FSV ORGANIC RESIN LENSES
INDEX 1.56 PREMIUM**



 ZhenJiang Pure OPTik Co., LTD

 +86 15952948899

 info@pureoptik.com

 www.pureoptik.com

Tag Description

 Premium	PREMIUM PRODUCT
 Basic	BASIC PRODUCT
 On Demand	ON DEMAND PRODUCT
 7 days	PRODUCT AVAILABLE IN 7 DAYS
	ON STOCK ITEM
	ON DEMAND ITEM
	RX ITEM
	ITEM AVAILABLE IN 7 DAYS
	SUITABLE FOR DAY DRIVING
	SUITABLE FOR NIGHT DRIVING
	NOT SUITABLE FOR NIGHT DRIVING

FSV Organic Resin Lens

INDEX 1,56 PREMIUM

HC / HMC Backside



GENERAL PRODUCT CHARACTERISTICS

INDEX	1.56
MATERIAL	KOC55/SR-56
DESIGN	Spheric/Aspherical
ABBE	35/41
DENSITY (gr/cm³)	1.27/1.15
COATING	HC / HTC/HMC/ SHMC Backside
COATING COLOR	Green/Blue/Revo/Anti-Glare
UV CUT (nm)	380/400/420++
FILTER UV CATEGORY	3 
PROGRESSIVE	Yes
BIOFOCAL	Yes
POLARIZED	Yes
TINTABLE	Yes
PHOTOCHROMIC	Matirial/Spin
SEMI-FINISHED	Yes
BAYER TEST	N/A

Pure Family

Pure

Poto



Your Optical Lenses Supply Chain Solution Provider

Available Range

[illegible]

AVAILABLE COLORS

●

