



**FINISHED
ORGANIC RESIN LENSES
TECHNICAL DATA**

**FSV ORGANIC RESIN LENSES
INDEX 1.57 PREMIUM**



 ZhenJiang Pure OPTik Co., LTD

 +86 15952948899

 info@pureoptik.com

 www.pureoptik.com

Tag Description

 Premium	PREMIUM PRODUCT
 Basic	BASIC PRODUCT
 On Demand	ON DEMAND PRODUCT
 7 days	PRODUCT AVAILABLE IN 7 DAYS
	ON STOCK ITEM
	ON DEMAND ITEM
	RX ITEM
	ITEM AVAILABLE IN 7 DAYS
	SUITABLE FOR DAY DRIVING
	SUITABLE FOR NIGHT DRIVING
	NOT SUITABLE FOR NIGHT DRIVING

FSV Organic Resin Lens

INDEX 1,57 Hi-impact PREMIUM



GENERAL PRODUCT CHARACTERISTICS	
INDEX	1.556 Hi-impact
MATERIAL	Hi-impact
DESIGN	Aspherical
ABBE	42
DENSITY (gr/cm³)	1.16
COATING	HMC /SHMC Backside
COATING COLOR	Green
UV CUT (nm)	400/420++
FILTER UV CATEGORY	3 
PROGRESSIVE	NO
BIOFOCAL	NO
POLARIZED	NO
TINTABLE	NO
PHOTOCHROMIC	NO
SEMI-FINISHED	NO
BAYER TEST	N/A

Pure Family

Pure

Poto



Your Optical Lenses Supply Chian Solution Provider

Available Range

Finished Single Vision 1.57 Hi-impact																																				
SPH	CYL(-,-)										SPH	CYL(+,-)																								
	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25		2.50	2.75	3.00	3.25	3.50	3.75	4.00	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00	
0.00																		0.00	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
0.25																		0.25																		
0.50																		0.50																		
0.75																		0.75																		
1.00																		1.00																		
1.25																		1.25																		
1.50																		1.50																		
1.75																		1.75																		
2.00																		2.00																		
2.25																		2.25																		
2.50																		2.50																		
2.75																		2.75																		
3.00																		3.00																		
3.25																		3.25																		
3.50																		3.50																		
3.75																		3.75																		
4.00																		4.00																		
4.25																		4.25																		
4.50																		4.50																		
4.75																		4.75																		
5.00																		5.00																		
5.25																		5.25																		
5.50																		5.50																		
5.75																		5.75																		
6.00																		6.00																		
6.25																																				
6.50																																				
6.75																																				
7.00																																				
7.25																																				
7.50																																				
7.75																																				
8.00																																				