



**FINISHED
ORGANIC RESIN LENSES
TECHNICAL DATA**

**FSV ORGANIC RESIN LENSES
INDEX 1.591 Polycarbonate PREMIUM**



 ZhenJiang Pure OPTik Co., LTD

 +86 15952948899

 info@pureoptik.com

 www.pureoptik.com

Tag Description

 Premium	PREMIUM PRODUCT
 Basic	BASIC PRODUCT
 On Demand	ON DEMAND PRODUCT
 7 days	PRODUCT AVAILABLE IN 7 DAYS
	ON STOCK ITEM
	ON DEMAND ITEM
	RX ITEM
	ITEM AVAILABLE IN 7 DAYS
	SUITABLE FOR DAY DRIVING
	SUITABLE FOR NIGHT DRIVING
	NOT SUITABLE FOR NIGHT DRIVING

FSV Organic Resin Lens

INDEX 1.59 Polycarbonate PREMIUM



GENERAL PRODUCT CHARACTERISTICS

INDEX	1.591
MATERIAL	Polycarbonate
DESIGN	Spheric
ABBE	32
DENSITY (gr/cm ³)	1.28
COATING	HC / HMC/SHMC Backside
COATING COLOR	Green/Blue
UV CUT (nm)	400/420++
FILTER UV CATEGORY	3
PROGRESSIVE	Yes
BIOFOCAL	NO
POLARIZED	NO
TINTABLE	HTC
PHOTOCHROMIC	Spin
SEMI-FINISHED	Yes
BAYER TEST	N/A

Pure Family

Pure

Poto



Your Optical Lenses Supply Chian Solution Provider

Available Range

Finished Single Vision 1.591 PC																																							
SPH	CYL(-,-)										SPH	CYL(+,-)																											
	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25		2.50	2.75	3.00	3.25	3.50	3.75	4.00	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00				
0.00	75mm										70mm										0.00	-	65mm										65mm						
0.25																					0.00	-																	
0.50																					0.00	-																	
0.75																					0.00	-																	
1.00																					0.00	-																	
1.25																					0.00	-																	
1.50																					0.00	-																	
1.75																					0.00	-																	
2.00																					0.00	-																	
2.25																					0.00	-																	
2.50																					0.00	-																	
2.75																					0.00	-																	
3.00	0.00	-																																					
3.25	70mm										65mm										3.25	-	65mm																
3.50																					3.00	-																	
3.75																					3.00	-																	
4.00																					3.00	-																	
4.25																					3.00	-																	
4.50																					3.00	-																	
4.75																					3.00	-																	
5.00																					3.00	-																	
5.25																					3.00	-																	
5.50																					3.00	-																	
5.75																					3.00	-																	
6.00																					3.00	-																	
6.25	70mm																				6.25	-																	
6.50																					6.00	-																	
6.75																					6.00	-																	
7.00																					6.00	-																	
7.25																					6.00	-																	
7.50																					6.00	-																	
7.75																					6.00	-																	
8.00																					6.00	-																	