



**FINISHED
ORGANIC RESIN LENSES
TECHNICAL DATA**

**FSV ORGANIC RESIN LENSES
INDEX 1.60 PREMIUM**



 ZhenJiang Pure OPTik Co., LTD

 +86 15952948899

 info@pureoptik.com

 www.pureoptik.com

Tag Description

 Premium	PREMIUM PRODUCT
 Basic	BASIC PRODUCT
 On Demand	ON DEMAND PRODUCT
 7 days	PRODUCT AVAILABLE IN 7 DAYS
	ON STOCK ITEM
	ON DEMAND ITEM
	RX ITEM
	ITEM AVAILABLE IN 7 DAYS
	SUITABLE FOR DAY DRIVING
	SUITABLE FOR NIGHT DRIVING
	NOT SUITABLE FOR NIGHT DRIVING

FSV Organic Resin Lens

INDEX 1,6 MR-8 PREMIUM



GENERAL PRODUCT CHARACTERISTICS

INDEX	1.60
MATERIAL	MR-8
DESIGN	spheric/Aspheric
ABBE	42
DENSITY (gr/cm ³)	1.30
COATING	HC / HTC/HMC/ SHMC Backside
COATING COLOR	Green/Blue
UV CUT (nm)	400/420++
FILTER UV CATEGORY	3 
PROGRESSIVE	Only SF
BIOFOCAL	Only Round Top
POLARIZED	Only SF
TINTABLE	Yes
PHOTOCHROMIC	Spin
SEMI-FINISHED	Yes
BAYER TEST	N/A

Pure Family

Pure

Poto



Your Optical Lenses Supply Chian Solution Provider

Available Range

Finished Single Vision 1.6 MR-8																													
SPH	CYL(-,-)																SPH	CYL(+,-)											
	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75		4.00	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00		
0.00	75mm									70mm								0.00	65mm										
0.25																		0.25											
0.50																		0.50											
0.75																		0.75											
1.00																		1.00											
1.25																		1.25											
1.50																		1.50											
1.75																		1.75											
2.00																		2.00											
2.25																		2.25											
2.50																		2.50											
2.75																		2.75											
3.00																		3.00											
3.25																		3.25											
3.50																		3.50											
3.75																		3.75											
4.00																		4.00											
4.25	4.25																												
4.50	4.50																												
4.75	4.75																												
5.00	5.00																												
5.25	5.25																												
5.50	5.50																												
5.75	5.75																												
6.00	6.00																												
6.25	70mm																												
6.50																													
6.75																													
7.00																													
7.25																													
7.50	70mm																												
7.75																													
8.00																													
8.25																													
8.50																													
8.75	70mm																												
9.00																													
9.25																													
9.50																													
9.75																													
10.00																													
10.25																													
10.50																													
10.75																													
11.00																													
11.25																													
11.50																													
11.75																													
12.00																													

AVAILABLE COLORS

- BROWN 469 
- GRAY 434 
- G15 343 