

Material GUIDE

ZhenJiang Pure Optik Co., LTD

Donggu Grou, Gucun Village, Situ Town, Danyang City, JIANGSU P.R., CHINA
Email: info@pureoptik.com Tel: +86 15952948899

Material List		
1	Product Name	1.499 Lens
	Material	CR-39
	Manufactory	PPG(USA); Acomon(Great LakeItaly)
	Refravtive Index (25℃)	n(d)=1.499; n(e)=1.501
	Abbe Value	58
	Specific Gravity (20℃)	1.32
2	Product Name	1.5 PHOTOCROMIC Lens
	Material	Acrylic resin
	Manufactory	Hgdbhn(Korea)
	Refravtive Index (25℃)	n(d)=1.503; n(e)=1.505
	Abbe Value	58
	Specific Gravity (20℃)	1.23
3	Product Name	Middle Index Lens
	Material	KOC 55
	Manufactory	KOC (Korea)
	Refravtive Index (25℃)	n(d)=1.542; n(e)=1.545
	Abbe Value	35
	Specific Gravity (20℃)	1.27
4	Product Name	PHOTOCROMIC Lens
	Material	SR-56
	Manufactory	Sartomer(USA)
	Refravtive Index (25℃)	n(d)=1.542; n(e)=1.545
	Abbe Value	41
	Specific Gravity (20℃)	1.15
5	Product Name	1.60 Lens
	Material	MR8
	Manufactory	Mitsui (Japan)
	Refravtive Index (25℃)	n(d)=1.594; n(e)=1.597
	Abbe Value	41
	Specific Gravity (20℃)	1.3
6	Product Name	1.67 Lens
	Material	MR7
	Manufactory	Mitsui (Japan)
	Refravtive Index (25℃)	n(d)=1.660; n(e)=1.665
	Abbe Value	31
	Specific Gravity (20℃)	1.35
7	Product Name	1.50 UV++Lens
	Material	Poly allyl carbonate ester
	Manufactory	Acomon(Italy)
	Refravtive Index (25℃)	n(d)=1.497; n(e)=1.498
	Abbe Value	58
	Specific Gravity (20℃)	1.28
8	Product Name	1.74 Lens
	Material	MR174
	Manufactory	Mitsui (Japan)
	Refravtive Index (25℃)	n(d)=1.732; n(e)=1.738
	Abbe Value	33
	Specific Gravity (20℃)	1.47