



PUREOPTIK

— Burn For Challenge —

**FINISHED
ORGANIC RESIN LENSES
TECHNICAL DATA**

**FSV ORGANIC RESIN LENSES
INDEX 1.5 PREMIUM**



 ZhenJiang Pure OPTik Co., LTD

 +86 15952948899

 info@pureoptik.com

 www.pureoptik.com

Tag Description



Premium

PREMIUM PRODUCT



Basic

BASIC PRODUCT



On Demand

ON DEMAND PRODUCT



7 days

PRODUCT AVAILABLE IN 7 DAYS



ON STOCK ITEM



ON DEMAND ITEM



RX ITEM



ITEM AVAILABLE IN 7 DAYS



SUITABLE FOR DAY DRIVING



SUITABLE FOR NIGHT DRIVING



NOT SUITABLE FOR NIGHT DRIVING

FSV Organic Resin Lens Index 1,5 CR-39 PREMIUM



GENERAL PRODUCT CHARACTERISTICS

INDEX	1,499
MATERIAL	CR-39
DESIGN	Spheric
ABBE	58
DENSITY (gr/cm ³)	1,32
COATING	HC / HTC/HMC/ SHMC Backside
COATING COLOR	Green/Blue/Revo/Anti-Glare
UV CUT (nm)	350/400/420++
FILTER UV CATEGORY	3 
PROGRESSIVE	Yes
BIOFOCAL	Yes
POLARIZED	Yes
TINTABLE	Yes
PHOTOCHROMIC	Spin
SEMI-FINISHED	Yes
BAYER TEST	N/A

Pure Family

Pure

Poto

Available Range

Finished Single Vision 1.50																																																																																				
SPH	CYL(-,-)																SPH	CYL(+,-)																																																																		
	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75		4.00	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00																																																	
0.00	65/70/75mm																0.00	50/55/60 65/70/75mm																55/60 65/70mm																																																		
0.25																																																																																				
0.50																																																																																				
0.75																																																																																				
1.00																																																																																				
1.25																																																																																				
1.50																																																																																				
1.75																																																																																				
2.00																																																																																				
2.25																																																																																				
2.50																																																																																				
2.75																																																																																				
3.00																																																																																				
3.25																																																																																				
3.50																																																																																				
3.75																																																																																				
4.00																																																																																				
4.25	65/70mm																4.25	50/55/60 65/70/75mm																																																																		
4.50																																																																																				
4.75																																																																																				
5.00																																																																																				
5.25																																																																																				
5.50																																																																																				
5.75																																																																																				
6.00																																																																																				
6.25	65/70mm																6.25	65mm																65mm																																																		
6.50																																																																																				
6.75																																																																																				
7.00																																																																																				
7.25																																																																																				
7.50																																																																																				
7.75																																																																																				
8.00																																																																																				

AVAILABLE COLORS

BROWN 469



GRAY 434



G15 343

