



PUREOPTIK

— Burn For Challenge —

**FINISHED
ORGANIC RESIN LENSES
TECHNICAL DATA**

**FSV ORGANIC RESIN LENSES
INDEX 1.56 PREMIUM**



 ZhenJiang Pure OPTik Co., LTD

 +86 15952948899

 info@pureoptik.com

 www.pureoptik.com

Tag Description

 Premium	PREMIUM PRODUCT
 Basic	BASIC PRODUCT
 On Demand	ON DEMAND PRODUCT
 7 days	PRODUCT AVAILABLE IN 7 DAYS
	ON STOCK ITEM
	ON DEMAND ITEM
	RX ITEM
	ITEM AVAILABLE IN 7 DAYS
	SUITABLE FOR DAY DRIVING
	SUITABLE FOR NIGHT DRIVING
	NOT SUITABLE FOR NIGHT DRIVING

FSV Organic Resin Lens INDEX 1,56 PREMIUM

HC / HMC Backside



GENERAL PRODUCT CHARACTERISTICS

INDEX	1.56
MATERIAL	KOC55/SR-56
DESIGN	Spheric/Aspherical
ABBE	35/41
DENSITY (gr/cm ³)	1.27/1.15
COATING	HC / HTC/HMC/ SHMC Backside
COATING COLOR	Green/Blue/Revo/Anti-Glare
UV CUT (nm)	380/400/420++
FILTER UV CATEGORY	3
PROGRESSIVE	Yes
BIOFOCAL	Yes
POLARIZED	Yes
TINTABLE	Yes
PHOTOCHROMIC	Matirial/Spin
SEMI-FINISHED	Yes
BAYER TEST	N/A

Pure Family

Pure

Poto

Available Range

Finished Single Vision 1.56																																																																																																			
SPH		CYL(-,-)																SPH		CYL(+,-)																																																																															
0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00																																																																		
0.00	72mm																0.00	55/60/65/70mm																0.00	55/60/65/70mm																																																																
0.25																	0.25																	70mm																	0.25	65mm																0.25	65mm																														
0.50																	0.50																																																																																		
0.75																	0.75																																																																																		
1.00																	1.00																																																																																		
1.25																	1.25																																																																																		
1.50																	1.50																																																																																		
1.75																	1.75																																																																																		
2.00																	2.00																																																																																		
2.25																	2.25																																																																																		
2.50																	2.50																																																																																		
2.75																	2.75																																																																																		
3.00	3.00																																																																																																		
3.25	3.25																																																																																																		
3.50	3.50																																																																																																		
3.75	3.75																																																																																																		
4.00	4.00																																																																																																		
4.25	4.25																																																																																																		
4.50	4.50																																																																																																		
4.75	4.75																																																																																																		
5.00	5.00																																																																																																		
5.25	5.25																																																																																																		
5.50	5.50																																																																																																		
5.75	5.75																																																																																																		
6.00	6.00																																																																																																		
6.25	6.25																																																																																																		
6.50	6.50																																																																																																		
6.75	6.75																																																																																																		
7.00	7.00																																																																																																		
7.25	7.25																																																																																																		
7.50	7.50																																																																																																		
7.75	7.75																																																																																																		
8.00	8.00																																																																																																		
8.25	70mm																8.25	65mm																8.25	65mm																																																																
8.50																	8.50																																																																																		
8.75																	8.75																																																																																		
9.00																	9.00																																																																																		
9.25	70mm																9.25	65mm																9.25	65mm																																																																
9.50																	9.50																																																																																		
9.75																	9.75																																																																																		
10.00																	10.00																																																																																		

AVAILABLE COLORS

