



PUREOPTIK

— Burn For Challenge —

**FINISHED
ORGANIC RESIN LENSES
TECHNICAL DATA**

**FSV ORGANIC RESIN LENSES
INDEX 1.57 PREMIUM**



 ZhenJiang Pure OPTik Co., LTD

 +86 15952948899

 info@pureoptik.com

 www.pureoptik.com

Tag Description



PREMIUM PRODUCT



BASIC PRODUCT



ON DEMAND PRODUCT



PRODUCT AVAILABLE IN 7 DAYS



ON STOCK ITEM



ON DEMAND ITEM



RX ITEM



ITEM AVAILABLE IN 7 DAYS



SUITABLE FOR DAY DRIVING



SUITABLE FOR NIGHT DRIVING



NOT SUITABLE FOR NIGHT DRIVING

FSV Organic Resin Lens INDEX 1,57 Hi-impact PREMIUM



GENERAL PRODUCT CHARACTERISTICS

INDEX	1.556 Hi-impact
MATERIAL	Hi-impact
DESIGN	Aspherical
ABBE	42
DENSITY (gr/cm³)	1.16
COATING	HMC /SHMC Backside
COATING COLOR	Green
UV CUT (nm)	400/420++
FILTER UV CATEGORY	3
PROGRESSIVE	NO
BIOFOCAL	NO
POLARIZED	NO
TINTABLE	NO
PHOTOCHROMIC	NO
SEMI-FINISHED	NO
BAYER TEST	N/A

Pure Family

Pure

Poto

Available Range

Finished Single Vision 1.57 Hi-impact																																																	
SPH	CYL(-,-)												SPH	CYL(+,-)																																			
	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75		3.00	3.25	3.50	3.75	4.00	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00														
0.00	75mm												0.00	65mm												65mm																							
0.25													0.25																																				
0.50													0.50																																				
0.75													0.75																																				
1.00													1.00																																				
1.25													1.25																																				
1.50													1.50																																				
1.75													1.75																																				
2.00													2.00																																				
2.25													2.25																																				
2.50													2.50																																				
2.75													2.75																																				
3.00													3.00																																				
3.25													3.25																																				
3.50													3.50																																				
3.75													3.75																																				
4.00													4.00																																				
4.25	70mm												4.25	65mm																																			
4.50													4.50																																				
4.75													4.75																																				
5.00													5.00																																				
5.25													5.25																																				
5.50													5.50																																				
5.75													5.75																																				
6.00													6.00																																				
6.25													70mm													6.25																							
6.50																																																	
6.75																																																	
7.00																																																	
7.25																																																	
7.50																																																	
7.75																																																	
8.00																																																	