



PUREOPTIK

— Burn For Challenge —

**FINISHED
ORGANIC RESIN LENSES
TECHNICAL DATA**

**FSV ORGANIC RESIN LENSES
INDEX 1.591 Polycarbonate PREMIUM**



 ZhenJiang Pure OPTik Co., LTD

 +86 15952948899

 info@pureoptik.com

 www.pureoptik.com

Tag Description



PREMIUM PRODUCT



BASIC PRODUCT



ON DEMAND PRODUCT



PRODUCT AVAILABLE IN 7 DAYS



ON STOCK ITEM



ON DEMAND ITEM



RX ITEM



ITEM AVAILABLE IN 7 DAYS



SUITABLE FOR DAY DRIVING



SUITABLE FOR NIGHT DRIVING



NOT SUITABLE FOR NIGHT DRIVING

FSV Organic Resin Lens INDEX 1.59 Polycarbonate PREMIUM



GENERAL PRODUCT CHARACTERISTICS

INDEX	1.591
MATERIAL	Polycarbonate
DESIGN	Spheric
ABBE	32
DENSITY (gr/cm ³)	1.28
COATING	HC / HMC/SHMC Backside
COATING COLOR	Green/Blue
UV CUT (nm)	400/420++
FILTER UV CATEGORY	3
PROGRESSIVE	Yes
BIOFOCAL	NO
POLARIZED	NO
TINTABLE	HTC
PHOTOCHROMIC	Spin
SEMI-FINISHED	Yes
BAYER TEST	N/A

Pure Family

Pure

Poto

Available Range

Finished Single Vision 1.591 PC																																					
SPH	CYL(-,-)												SPH	CYL(+,-)																							
	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75		3.00	3.25	3.50	3.75	4.00	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00		
0.00	75mm												0.00	65mm												65mm											
0.25																																					
0.50																																					
0.75																																					
1.00																																					
1.25																																					
1.50																																					
1.75																																					
2.00																																					
2.25																																					
2.50																																					
2.75																																					
3.00	70mm												3.00	65mm												65mm											
3.25																																					
3.50																																					
3.75																																					
4.00																																					
4.25																																					
4.50																																					
4.75																																					
5.00																																					
5.25																																					
5.50																																					
5.75																																					
6.00	70mm												6.00	65mm												65mm											
6.25																																					
6.50																																					
6.75																																					
7.00																																					
7.25																																					
7.50																																					
7.75																																					
8.00	70mm													65mm												65mm											