



PUREOPTIK

— Burn For Challenge —

**FINISHED
ORGANIC RESIN LENSES
TECHNICAL DATA**

**FSV ORGANIC RESIN LENSES
INDEX 1.60 PREMIUM**



 ZhenJiang Pure OPTik Co., LTD

 +86 15952948899

 info@pureoptik.com

 www.pureoptik.com

Tag Description



PREMIUM PRODUCT



BASIC PRODUCT



ON DEMAND PRODUCT



PRODUCT AVAILABLE IN 7 DAYS



ON STOCK ITEM



ON DEMAND ITEM



RX ITEM



ITEM AVAILABLE IN 7 DAYS



SUITABLE FOR DAY DRIVING



SUITABLE FOR NIGHT DRIVING



NOT SUITABLE FOR NIGHT DRIVING

FSV Organic Resin Lens INDEX 1,6 MR-8 PREMIUM



GENERAL PRODUCT CHARACTERISTICS

INDEX	1.60
MATERIAL	MR-8
DESIGN	spheric/Aspheric
ABBE	42
DENSITY (gr/cm ³)	1.30
COATING	HC / HTC/HMC/ SHMC Backside
COATING COLOR	Green/Blue
UV CUT (nm)	400/420++
FILTER UV CATEGORY	3
PROGRESSIVE	Only SF
BIOFOCAL	Only Round Top
POLARIZED	Only SF
TINTABLE	Yes
PHOTOCHROMIC	Spin
SEMI-FINISHED	Yes
BAYER TEST	N/A

Pure Family

Pure

Poto

Available Range

Finished Single Vision 1.6 MR-8																																																								
CYL(-,-)														CYL(+,-)																																										
SPH	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00	SPH	0.00	0.25	0.50	0.75	1.00	1.25	1.50	1.75	2.00																													
0.00	75mm														70mm														0.00	-	-	-	-	-	-	-	-	-	-	-	-	-	-													
0.25																													0.00	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	0.25	-	-	-	-	-	-	-	-	-	-	-
0.50																													0.50	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	0.50	-	-	-	-	-	-	-	-	-	-	-
0.75																													0.75	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	0.75	-	-	-	-	-	-	-	-	-	-	-
1.00																													1.00	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	1.00	-	-	-	-	-	-	-	-	-	-	-
1.25																													1.25	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	1.25	-	-	-	-	-	-	-	-	-	-	-
1.50																													1.50	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	1.50	-	-	-	-	-	-	-	-	-	-	-
1.75																													1.75	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	1.75	-	-	-	-	-	-	-	-	-	-	-
2.00																													2.00	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	2.00	-	-	-	-	-	-	-	-	-	-	-
2.25																													2.25	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	2.25	-	-	-	-	-	-	-	-	-	-	-
2.50																													2.50	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	2.50	-	-	-	-	-	-	-	-	-	-	-
2.75																													2.75	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	2.75	-	-	-	-	-	-	-	-	-	-	-
3.00																													3.00	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	3.00	-	-	-	-	-	-	-	-	-	-	-
3.25																													3.25	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	3.25	-	-	-	-	-	-	-	-	-	-	-
3.50																													3.50	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	3.50	-	-	-	-	-	-	-	-	-	-	-
3.75																													3.75	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	3.75	-	-	-	-	-	-	-	-	-	-	-
4.00	4.00	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	4.00	-	-	-	-	-	-	-	-	-	-	-																												
4.25	4.25	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	4.25	-	-	-	-	-	-	-	-	-	-	-																												
4.50	4.50	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	4.50	-	-	-	-	-	-	-	-	-	-	-																												
4.75	4.75	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	4.75	-	-	-	-	-	-	-	-	-	-	-																												
5.00	5.00	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	5.00	-	-	-	-	-	-	-	-	-	-	-																												
5.25	5.25	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	5.25	-	-	-	-	-	-	-	-	-	-	-																												
5.50	5.50	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	5.50	-	-	-	-	-	-	-	-	-	-	-																												
5.75	5.75	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	5.75	-	-	-	-	-	-	-	-	-	-	-																												
6.00	6.00	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	6.00	-	-	-	-	-	-	-	-	-	-	-																												
6.25	70mm																																																							
6.50																																																								
6.75																																																								
7.00																																																								
7.25																																																								
7.50																																																								
7.75																																																								
8.00																																																								
8.25	70mm																																																							
8.50																																																								
8.75																																																								
9.00																																																								
9.25																																																								
9.50																																																								
9.75																																																								
10.00																																																								
10.25	70mm																																																							
10.50																																																								
10.75																																																								
11.00																																																								
11.25																																																								
11.50																																																								
11.75																																																								
12.00																																																								

AVAILABLE COLORS

BROWN 469



GRAY 434



G15 343

